



LABORATORY VERTICAL AUTOCLAVE BAVT-2901

LABORATORY VERTICAL AUTOCLAVE BAVT-2901

VERTICAL AUTOCLAVE



Three times pre-vacuum
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Safety interlocking device
LED touch screen display
Standard test interface
Self-expanding seal
With drying function
With stainless steel bucket or basket
With automatic protection function/over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning.
Automatic water

SPECIFICATIONS

Model	BAVT-2901
Capacity	35 L
Power	3.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Drying time setting range	0 ~ 240 min
Chamber volume	Ø 350x400
Bucket dimension	Ø 330x320
Basket dimension	Ø 320x350
Outer dimension	698x498x940
Packing size	790x610x1130
Net weight	115 Kg
Gross Weight	130 Kg
Alt Name	Vertical Steam Sterilizer

FEATURES

Three times pre-vacuum
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Safety interlocking device
LED touch screen display
Standard test interface
Self-expanding seal
With drying function
With stainless steel bucket or basket
With automatic protection function/over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning.
Automatic water

External steam Generator

Buzzer reminder after sterilization, automatic stop.

Automatic discharge of cold air, automatic steam exhaust after sterilization

Internal printer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com