



LABORATORY VERTICAL AUTOCLAVE BAVT-3004

LABORATORY VERTICAL AUTOCLAVE BAVT-3004

VERTICAL AUTOCLAVE



Hand wheel translation and quick opening structure
High quality stainless steel material
With drying function
Microcomputer automatic control, arbitrary setting sterilization parameters
Safety interlocking device
Digital display and Touch buttons
Standard test interface
Self-expanding seal
With automatic protection function over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning
With stainless steel bucket or basket
Buzzer reminder after sterilization, automatic stop

SPECIFICATIONS

Model	BAVT-3004
Capacity	100 L
Power	4.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Drying time setting range	0 ~ 240 min
Chamber volume	Ø 450x650
Bucket dimension	Ø 420x540
Basket dimension	Ø 410x300x2
Outer dimension	788x588x1255
Packing size	870x650x1380
Net weight	129 Kg
Gross Weight	156 Kg
Alt Name	Vertical Steam Sterilizer

FEATURES

Hand wheel translation and quick opening structure
High quality stainless steel material
With drying function
Microcomputer automatic control, arbitrary setting sterilization parameters
Safety interlocking device
Digital display and Touch buttons
Standard test interface
Self-expanding seal
With automatic protection function over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning
With stainless steel bucket or basket
Buzzer reminder after sterilization, automatic stop

Automatic discharge of cold air, automatic steam exhaust after sterilization, Internal steam cycle



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com