



LABORATORY VERTICAL AUTOCLAVE BAVT-3201

LABORATORY VERTICAL AUTOCLAVE BAVT-3201

VERTICAL AUTOCLAVE



Hand wheel translation and quick opening structure
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Safety interlocking device
Digital display and Touch buttons
Standard test interface
Self-expanding seal
With automatic protection function: over temperature protection; over pressure self discharge protection; low water level protection; anti dry burning
With stainless steel bucket or basket
Buzzer reminder after sterilization, automatic stop
Automatic discharge of cold air, automatic steam exhaust after sterilization

SPECIFICATIONS

Model	BAVT-3201
Capacity	35 L
Power	3.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 350x400
Bucket dimension	Ø 330x320
Basket dimension	Ø 320x350
Outer dimension	480x480x1000
Packing size	580x630x1130
Net weight	58 Kg
Gross Weight	79 Kg
Alt Name	Vertical Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com