



LABORATORY VERTICAL AUTOCLAVE BAVT-3202

LABORATORY VERTICAL AUTOCLAVE BAVT-3202

VERTICAL AUTOCLAVE



Hand wheel translation and quick opening structure
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Safety interlocking device
Digital display and Touch buttons
Standard test interface
Self-expanding seal
With automatic protection function: over temperature protection; over pressure self discharge protection; low water level protection; anti dry burning
With stainless steel bucket or basket
Buzzer reminder after sterilization, automatic stop
Automatic discharge of cold air, automatic steam exhaust after sterilization

SPECIFICATIONS

Model	BAVT-3202
Capacity	50 L
Power	3.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 350x525
Bucket dimension	Ø 330x480
Basket dimension	Ø 320x240x2
Outer dimension	480x480x1120
Packing size	580x630x1230
Net weight	80 Kg
Gross Weight	106 Kg
Alt Name	Vertical Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com