



# LABORATORY VERTICAL AUTOCLAVE BAVT-3204

# LABORATORY VERTICAL AUTOCLAVE BAVT-3204

## VERTICAL AUTOCLAVE



Hand wheel translation and quick opening structure  
High quality stainless steel material  
Microcomputer automatic control, arbitrary setting sterilization parameters  
Safety interlocking device  
Digital display and Touch buttons  
Standard test interface  
Self-expanding seal  
With automatic protection function: over temperature protection; over pressure self discharge protection; low water level protection; anti dry burning  
With stainless steel bucket or basket  
Buzzer reminder after sterilization, automatic stop  
Automatic discharge of cold air, automatic steam exhaust after sterilization

## SPECIFICATIONS

Model	BAVT-3204
Capacity	100 L
Power	4.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 450x650
Bucket dimension	Ø 430x550
Basket dimension	Ø 410x300x2
Outer dimension	555x555x1250
Packing size	690x650x1410
Net weight	106 Kg
Gross Weight	133 Kg
Alt Name	Vertical Steam Sterilizer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)