



LABORATORY VERTICAL AUTOCLAVE BAVT-3301

LABORATORY VERTICAL AUTOCLAVE BAVT-3301

VERTICAL AUTOCLAVE



Bolt structure
High quality stainless steel material
Microcomputer automatic control, arbitrary set of sterilization parameters
Digital display and touch buttons
Standard test interface
The container with automatic protection function: over temperature protection;
over pressure protection self discharge; low water level protection, anti dry burning
With stainless steel bucket or basket
Buzzer reminder after sterilization, automatic stop.
Automatic cooling air, sterilization end automatic exhaust steam

SPECIFICATIONS

Model	BAVT-3301
Capacity	35 L
Power	2.5 kW
Voltage	AC 220V, 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 300x500
Bucket dimension	Ø 280x400
Basket dimension	/
Outer dimension	400x400x860
Packing size	550x500x1030
Net weight	41 Kg
Gross Weight	58 Kg
Alt Name	Vertical Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com