



# LABORATORY VERTICAL AUTOCLAVE BAVT-3302

# LABORATORY VERTICAL AUTOCLAVE BAVT-3302

## VERTICAL AUTOCLAVE



Bolt structure  
High quality stainless steel material  
Microcomputer automatic control, arbitrary set of sterilization parameters  
Digital display and touch buttons  
Standard test interface  
The container with automatic protection function: over temperature protection;  
over pressure protection self discharge; low water level protection, anti dry burning  
With stainless steel bucket or basket  
Buzzer reminder after sterilization, automatic stop.  
Automatic cooling air, sterilization end automatic exhaust steam

## SPECIFICATIONS

Model	BAVT-3302
Capacity	50 L
Power	3.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 350x525
Bucket dimension	Ø 330x460
Basket dimension	Ø 320x240x2
Outer dimension	480x480x960
Packing size	580x630x1100
Net weight	53 Kg
Gross Weight	75 Kg
Alt Name	Vertical Steam Sterilizer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)