



LABORATORY VERTICAL AUTOCLAVE BAVT-3304

LABORATORY VERTICAL AUTOCLAVE BAVT-3304

VERTICAL AUTOCLAVE



Bolt structure
High quality stainless steel material
Microcomputer automatic control, arbitrary set of sterilization parameters
Digital display and touch buttons
Standard test interface
The container with automatic protection function: over temperature protection;
over pressure protection self discharge; low water level protection, anti dry burning
With stainless steel bucket or basket
Buzzer reminder after sterilization, automatic stop.
Automatic cooling air, sterilization end automatic exhaust steam

SPECIFICATIONS

Model	BAVT-3304
Capacity	100 L
Power	4.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 450x650
Bucket dimension	Ø 420x540
Basket dimension	Ø 410x300x2
Outer dimension	555x555x1080
Packing size	690x650x1230
Net weight	78 Kg
Gross Weight	103 Kg
Alt Name	Vertical Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com