



# LABORATORY VERTICAL AUTOCLAVE BAVT-3305

# LABORATORY VERTICAL AUTOCLAVE BAVT-3305

## VERTICAL AUTOCLAVE



Bolt structure  
High quality stainless steel material  
Microcomputer automatic control, arbitrary set of sterilization parameters  
Digital display and touch buttons  
Standard test interface  
The container with automatic protection function: over temperature protection;  
over pressure protection self discharge; low water level protection, anti dry burning  
With stainless steel bucket or basket  
Buzzer reminder after sterilization, automatic stop.  
Automatic cooling air, sterilization end automatic exhaust steam

## SPECIFICATIONS

Model	BAVT-3305
Capacity	120 L
Power	4.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 450x760
Bucket dimension	Ø 420x660
Basket dimension	Ø 410x350x2
Outer dimension	555x599x1200
Packing size	690x650x1330
Net weight	82 Kg
Gross Weight	108 Kg
Alt Name	Vertical Steam Sterilizer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)