



## PRECISION BALANCE BBAL-209

# PRECISION BALANCE BBAL-209

PRECISION BALANCE



Fashional appearance  
 Large size stainless steel pan (For XY-1BF)  
 LCD (White back light with black font)  
 AC and DC power supply (AA\*3)  
 Tare function / Counting / Unit conversion (g/ct/oz...)  
 Min weighing set  
 Overload alarm / Level indicator  
 Option: Dual display / Interface / Printer / Dust cover  
 Calibration weight

## SPECIFICATIONS

Model	BBAL-209
Capacity	10000g
Readability	0.1g
Repeatability	±0.2g
Linearity	±0.3g
Pan Size	160x160mm
Stable Time	≤2S
Alt Name	Precision Balance

## ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g
2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g



1



1




1



1



2



2



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
 Email: info@biolabscientific.com | Website: www.biolabscientific.com