



## PRECISION BALANCE BBAL-217

# PRECISION BALANCE BBAL-217

PRECISION BALANCE



LCD / AC and DC power supply (AA\*3)  
 Tare function / Counting / Unit conversion (g/ct/oz...)  
 Min weighing set  
 Overload alarm / Level indicator  
 Option: Dual display / Interface / Printer / Windproof / Dust cover  
 Note: Capacity more than 1 kg or accuracy less than 0.01 g is not equipped with calibration weight.

## SPECIFICATIONS

Model	BBAL-217
Capacity	310g
Readability	0.01g
Repeatability	±0.02g
Linearity	±0.03g
Pan Size	Φ130mm
Stable Time	≤2S
Alt Name	Precision Balance

## ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g
2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g



1



1



1



1




2




2

No	Name	Size
1	Plastic windshield	18*17*9.5cm
2	Collapsible plastic windshield	18*18.5*9cm
3	Dust cover	19*27*6.5cm



1



2



3



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)