



## PRECISION BALANCE BBAL-220

# PRECISION BALANCE BBAL-220

PRECISION BALANCE



LCD/AC and DC power supply (AA\*3)  
 Tare function / Counting / Unit conversion (g/ct/oz...)  
 Min weighing set  
 Overload alarm / Level indicator  
 Option: Dual display / Interface / Printer / Dust cover

## SPECIFICATIONS

Model	BBAL-220
Capacity	110g
Readability	0.001g
Min Weighing	0.001g-0.004g
Repeatability	±0.002g
Linearity	±0.003g
Stable Time	≤3s
Pan Size	Φ90mm
Appearance Size	270*200*280mm
Draft Shield Size	170*165*180mm (Inner diameter)
Calibration	External calibration
Cal. Weight	100g
Packing Size	360*265*350mm / 560*380*730mm (4 units per box)
Gross Weight	4kg / 16kg (4 units per box)
Alt Name	Precision Balance

## ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g
2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g



1



1



1



1



2



2



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)