



## PRECISION BALANCE BBAL-224

# PRECISION BALANCE BBAL-224

PRECISION BALANCE



Aluminum shell  
 Unit conversion (g/ct/oz...)  
 LED (Green) / LCD  
 Tare function / Counting  
 Overload alarm / Level indicator  
 Option: Printer

## SPECIFICATIONS

|               |                   |
|---------------|-------------------|
| Model         | BBAL-224          |
| Capacity      | 4000g             |
| Readability   | 0.01g             |
| Repeatability | ±0.02g            |
| Linearity     | ±0.03g            |
| Pan Size      | 180*180mm         |
| Stable Time   | ≤3s               |
| Interface     | RS232 (Standard)  |
| Dimensions    | 345*205*105mm     |
| Alt Name      | Precision Balance |

## ACCESSORIES FOR PURCHASE

| No | Name       | Description   | Information  |
|----|------------|---|--|
| 1  | Weight set | Recommend corresponding accuracy calibration (1mg-20kg) | M1 level Suitable for 0.1g-0.01g<br>F2 level Suitable for 0.01g-0.001g<br>F1 level Suitable for 0.001g-0.0001g<br>E2 level Suitable for over 0.0001g |
| 2  | Weight set | Lock type series weights (5kg-20kg)                     | M1 level Suitable for 0.1g-0.01g<br>F2 level Suitable for 0.01g-0.001g<br>F1 level Suitable for 0.001g-0.0001g                                       |



1



1



1



1



2



2

## FEATURES

Aluminum shell  
 Unit conversion (g/ct/oz...)  
 LED (Green) / LCD

Tare function / Counting  
Overload alarm / Level indicator  
Option: Printer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)