



## PRECISION BALANCE BBAL-225

# PRECISION BALANCE BBAL-225

PRECISION BALANCE



Aluminum shell  
 Unit conversion (g/ct/oz...)  
 LED (Green) / LCD  
 Tare function / Counting  
 Overload alarm / Level indicator  
 Option: Printer

## SPECIFICATIONS

Model	BBAL-225
Capacity	5000g
Readability	0.01g
Repeatability	±0.02g
Linearity	±0.03g
Pan Size	180*180mm
Stable Time	≤3s
Interface	RS232 (Standard)
Dimensions	345*205*105mm
Alt Name	Precision Balance

## ACCESSORIES FOR PURCHASE

No	Name	Description	Information
1	Weight set	Recommend corresponding accuracy calibration (1mg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g E2 level Suitable for over 0.0001g
2	Weight set	Lock type series weights (5kg-20kg)	M1 level Suitable for 0.1g-0.01g F2 level Suitable for 0.01g-0.001g F1 level Suitable for 0.001g-0.0001g



1



1



1



1



2



2

## FEATURES

Aluminum shell  
 Unit conversion (g/ct/oz...)  
 LED (Green) / LCD

Tare function / Counting  
Overload alarm / Level indicator  
Option: Printer



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)