



## WEIGHING SCALE BBAL-402

# WEIGHING SCALE BBAL-402

INDUSTRIAL SCALE



Large range and high precision  
Fashionable appearance  
Large size stainless steel pan  
LCD (White back light with black font)  
AC and DC power supply  
Tare function / Counting /  
Unit conversion (g/ct/oz), Option lb  
Overload alarm / Level indicator  
Option: Interface / Printer

## SPECIFICATIONS

Model	BBAL-402
Readability	0.1g
Capacity	5.1kg
Repeatability	±0.2g
Linearity	±0.3g
Pan Size	315x230mm
Stable Time	≤2S
Alt Name	Industrial scale



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)