



WEIGHING SCALE BBAL-404

WEIGHING SCALE BBAL-404

INDUSTRIAL SCALE



Large range and high precision
Fashionable appearance
Large size stainless steel pan
LCD (White back light with black font)
AC and DC power supply
Tare function / Counting /
Unit conversion (g/ct/oz), Option lb
Overload alarm / Level indicator
Option: Interface / Printer

SPECIFICATIONS

Model	BBAL-404
Readability	0.1g
Capacity	16kg
Repeatability	±0.2g
Linearity	±0.3g
Pan Size	315x230mm
Stable Time	≤2S
Alt Name	Industrial scale



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com