



WEIGHING SCALE BBAL-414

WEIGHING SCALE BBAL-414


INDUSTRIAL SCALE



Large range and high precision
Spray steel scale frame
Large size stainless steel table
LED (Red)
AC and DC power supply
Tare function / Counting / Unit conversion (kg/lb)
Overload alarm / Undervoltage alarm / Level indicator
Option: Interface / Printer

SPECIFICATIONS

Model	BBAL-414
Capacity	310kg
Readability	10g
Repeatability	±20g
Linearity	±30g
Pan Size	400x500mm
Stable Time	≤2S
Alt Name	Industrial scale



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com