



## WEIGHING SCALE BBAL-418

# WEIGHING SCALE BBAL-418

INDUSTRIAL SCALE



Large range and high precision  
Spray steel scale frame  
Large size stainless steel table  
LED (Red)  
AC and DC power supply  
Tare function / Counting / Unit conversion (kg/lb)  
Overload alarm / Undervoltage alarm / Level indicator  
Option: Interface / Printer

## SPECIFICATIONS

Model	BBAL-418
Capacity	110kg
Readability	1g
Repeatability	±2g
Linearity	±3g
Pan Size	400x500mm
Stable Time	≤2S
Alt Name	Industrial scale



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)