



BLOOD BANK REFRIGERATED CENTRIFUGE BBRC-202

BLOOD BANK REFRIGERATED CENTRIFUGE BBRC-202

BLOOD BANK REFRIGERATED CENTRIFUGE



- 1.Floor type with high quality steel structure,safe and reliable.
- 2.Flexible axle driven system which drive the rotor directly, smooth in operation, low noise.
- 3.Brushless DC frequency motor with simpler construction,more reliable performance,longer life and quietly running
- 2.4.Imported compressors,fluorine-free, double cycle cooling, cold and hot alternating easily, environment pollution-free and precise in temperature control.
- 4.Microprocessor control,digital display in the speed, time, RCF in operation, speed rising and reducing quick, operate simply.
- 5.10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating.
- 6.Automatically electric lid lock, super speed, over temperature protection and

SPECIFICATIONS

Model	BBRC-202
Max speed	6000 rpm
Max volume	6x1000ml (6 bags)
Timer	1~9h59min
Dimension(HxDxW)	800x740x930mm
Speed accuracy	±20rpm
Temperature accuracy	-20~40°C
Max RCF	6600xg
Noise	≤60dBA
Net weight	260KG
Power supply	AC 220V 50Hz 20A
Temperature range	±1°C
Package	wooden box
Alt Name	Blood Bank Refrigerated Centrifuge

ITEMS INCLUDED

No	Name	Max Speed (rpm)	Max Volume (ml)	Max RCF (g)
1	Swing Rotor	4200	6x1000ml(6 bags)	5100



1

FEATURES

- 1.Floor type with high quality steel structure,safe and reliable.
- 2.Flexible axle driven system which drive the rotor directly, smooth in operation, low noise.
- 3.Brushless DC frequency motor with simpler construction,more reliable performance,longer life and quietly running
- 2.4.Imported compressors,fluorine-free, double cycle cooling, cold and hot alternating easily, environment pollution-free and precise in temperature control.
- 4.Microprocessor control,digital display in the speed, time, RCF in operation, speed rising and reducing quick, operate simply.

- 5.10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating.
6. Automatically electric lid lock, super speed, over temperature protection and imbalance protection.
7. Three layer protection of steel structure, safe and reliable.
8. It is consistent with CE security certification and ISO9001 quality certification.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com