



BLOOD BANK REFRIGERATED CENTRIFUGE BBRC-205

BLOOD BANK REFRIGERATED CENTRIFUGE BBRC-205

BLOOD BANK REFRIGERATED CENTRIFUGE



1. Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running
2. Flexible axle driven system which drive the rotor directly, smooth in operation, low noise.
3. Imported compressors, fluorine-free, double cycle cooling, cold and hot alternating easily, environment pollution-free and precise in temperature control.
4. Microprocessor control, digital display in the speed, time, RCF in operation, speed rising and reducing quick, operate simply.
5. 10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating.
6. Automatically electric lid lock, super speed, over temperature protection and imbalance protection. The centrifuge body is made of high quality steel, safe an

SPECIFICATIONS

Model	BBRC-205
Max speed	7000rpm
Max volume	6x2400ml(12 bags)
Timer	1~9h59min
Dimension(HxDxW)	940x890x1000mm
Speed accuracy	±20rpm
Temperature accuracy	±1° C
Max RCF	11900xg
Noise	≤60dBA
Net weight	560KG
Power supply	AC 380V 50HZ 35A
Temperature range	-20~40° C
Package	wooden box
Alt Name	Blood Bank Refrigerated Centrifuge

ITEMS INCLUDED

No	Name	Max Speed (rpm)	Max Volume (ml)	Max RCF (g)
1	Swing Rotor	4000	6x2400ml(12 bags)	5530



1

FEATURES

1. Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running
2. Flexible axle driven system which drive the rotor directly, smooth in operation, low noise.
3. Imported compressors, fluorine-free, double cycle cooling, cold and hot alternating easily, environment pollution-free and precise in temperature control.
4. Microprocessor control, digital display in the speed, time, RCF in operation, speed rising and reducing quick, operate simply.
5. 10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating.

6. Automatically electric lid lock, super speed, over temperature protection and imbalance protection. The centrifuge body is made of high quality steel, safe and reliable.
7. Three layer protection of steel structure, safe and reliable.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com