



PRODUCT CATALOG



GEL INSTRUMENT BBRC-302

GEL INSTRUMENT BBRC-302

GEL CARD INCUBATOR



1. Keep 24-Gel Cards in temperature of 37 degree.
2. 2 independent incubation zones with a capacity for 12 Gel cards
3. Digital display to set up time and keep the temperature
4. Audible programmable alarm to notify end of incubation period.

SPECIFICATIONS

Model	BBRC-302
Ambient Temperature	10~37°C
Temperature uniformity	≤±0.5°C
Max. Capacity	24 Blood cards
Size (W x D x H mm)	390 x 280 x 140 mm
Voltage(V/Hz)	AC 220V/110V 50HZ/60HZ
Incubation Temperature	37°C ±0.1°C
Temperature Accuracy	<0.5°C
Rated Power	250W
Net Weight	6.2 KG
Time Range	0~99min59sec
Certificates	CE, ISO & Calibration report are available
Alt Name	Gel Card Incubator

APPLICATIONS

Digital incubator for the Gel cards. Blood card/ID card incubator holds 2 independent incubation zones with a capacity for 12 Gel cards. It can be suitable for many types of gel cards from different manufacturers.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com