



## SHAKING WATER BATH BBSW-103

# SHAKING WATER BATH BBSW-103

## SHAKING WATER BATH



Microcomputer temperature and frequency controller, LCD display screen;  
Stainless steel chamber;  
With power outage recovery and parameter memory function;  
Universal spring fixture;  
Equipped with timing function.

## SPECIFICATIONS

Model	BBSW-103
Type	Water Bath
Temp. Range [°C]	RT+5~100
Oscillation Mode	Orbital or Reciprocating
Rotation Rate [rpm]	0~220 (Orbital) 0~300 (Reciprocating)
Vibration Amplitude [mm]	20 (Orbital) 30 (Reciprocating)
Power [W]	1600
Timer Range [min/h]	0~9999
Internal Dimensions WxDxH [mm]	480x390x170
Housing Dimensions WxDxH [mm]	705x515x470
Tray Size WxDxH [mm]	380x300x170
Full Load Configuration [mlxpcs]	500x6 or 1000x3
Net Weight/Gross weight [kg]	50/68
Package Size [mm]	810x600x600
Alt Name	Shaking Water Bath



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)