



SHAKING WATER BATH BBSW-3004

SHAKING WATER BATH BBSW-3004

SHAKING WATER BATH

LED Microprocessor Controller (with timing function)

Widely applicable for laboratory researches on bacteria cultivation, fermentation, hybridization, chemical and biochemical reaction, enzymes and tissues research, which have a high requirement on precision of shaking speed and temperature.



SPECIFICATIONS

Model	BBSW-3004
Temperature Range	RT+5~100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking Speed Range	30~180rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1650W
Interior Dimension (WxDxH) mm	618x310x250
Exterior Dimension (WxDxH) mm	823x350x355
Packing Dimension	48*58*94
CBM	0.26
Net Weight	29.00
Gross weight	45
Alt Name	Shaking Water Bath





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com