



SHAKING WATER BATH BBSW-3005

SHAKING WATER BATH BBSW-3005

SHAKING WATER BATH

LED Microprocessor Controller (with timing function)

Widely applicable for laboratory researches on bacteria cultivation, fermentation, hybridization, chemical and biochemical reaction, enzymes and tissues research, which have a high requirement on precision of shaking speed and temperature.



SPECIFICATIONS

Model	BBSW-3005
Temperature Range	10~100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking Speed Range	30~180rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1500W
Interior Dimension (WxDxH) mm	440x300x250
Exterior Dimension (WxDxH) mm	710x410x710
Packing Dimension	84*54*89
CBM	0.40
Net Weight	60.00
Gross weight	83
Alt Name	Shaking Water Bath





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com