



WATER BATH BBWA-3101

WATER BATH BBWA-3101

WATER BATH



Microprocessor controller (with timing function)

- A stamping molding stainless steel tank, easy to clean.
- LCD screen, multiple data display with timing function, easy to operate.
- Stainless steel shelves cover heater and sensor to avoid damage during using.
- Once-forming stainless steel lid.
- Cut off heater automatically in case of lack of water, meanwhile visible and audible alarms to remind users in time.
- Independent temperature-limiting alarm system.
- temperature error alarm.
- Test tube holder can be placed. (Option)

SPECIFICATIONS

Model	BBWA-3101
Electrical Requirement	AC220 50Hz
Power Consumption	400W
Temperature Range	RT+5~100°C
Temperature Stability	±0.3°C
Temp Alarm	±2°C
Interior Dimension (WxDxH) mm	130x280x150
Exterior Dimension (WxDxH) mm	396x250x260
Timing Range	1~5999min
Chamber Volume	2 holes Φ 112mm
Packing Dimension	52*29*30
CBM	0.05
Net Weight	0.05
Gross weight	8
Alt Name	Water Bath



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com