



## WATER BATH BBWA-3105

# WATER BATH BBWA-3105

WATER BATH



Microprocessor controller (with timing function)

- A stamping molding stainless steel tank, easy to clean.
- LCD screen, multiple data display with timing function, easy to operate.
- Stainless steel shelves cover heater and sensor to avoid damage during using.
- Once-forming stainless steel lid.
- Cut off heater automatically in case of lack of water, meanwhile visible and audible alarms to remind users in time.
- Independent temperature-limiting alarm system.
- temperature error alarm.
- Test tube holder can be placed. (Option)

## SPECIFICATIONS

Model	BBWA-3105
Electrical Requirement	AC220 50Hz
Power Consumption	1600W
Temperature Range	RT+5~100°C
Temperature Stability	±0.3°C
Temp Alarm	±2°C
Interior Dimension (WxDxH) mm	280x220x150
Exterior Dimension (WxDxH) mm	526x395x260
Timing Range	1~5999min
Chamber Volume	2 holes + 4 holes
Packing Dimension	60*47*33
CBM	0.09
Net Weight	13.00
Gross weight	15
Alt Name	Water Bath



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)