



WATER BATH BBWA-3106

WATER BATH BBWA-3106

WATER BATH



Microprocessor controller (with timing function)

- A stamping molding stainless steel tank, easy to clean.
- LCD screen, multiple data display with timing function, easy to operate.
- Stainless steel shelves cover heater and sensor to avoid damage during using.
- Once-forming stainless steel lid.
- Cut off heater automatically in case of lack of water, meanwhile visible and audible alarms to remind users in time.
- Independent temperature-limiting alarm system.
- temperature error alarm.
- Test tube holder can be placed. (Option)

SPECIFICATIONS

| | |
|-------------------------------|-------------|
| Model | BBWA-3106 |
| Electrical Requirement | AC220 50Hz |
| Power Consumption | 800W |
| Temperature Range | RT+5~100°C |
| Temperature Stability | ±0.2°C |
| Temp Alarm | 0.1°C |
| Interior Dimension (WxDxH) mm | 300x240x200 |
| Exterior Dimension (WxDxH) mm | 480x300x480 |
| Timing Range | 1~5999min |
| Chamber Volume | 11L |
| Packing Dimension | 57*40*47 |
| CBM | 0.11 |
| Net Weight | 12.00 |
| Gross weight | 15 |
| Alt Name | Water Bath |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com