



## WATER BATH BBWA-3201

# WATER BATH BBWA-3201

CIRCULATING BATH, WATER BATH



- Microprocessor controller (with timing function)
- Microprocessor control with timing function.
  - Digital display.
  - Audible and visible alarm for over temperature.

## SPECIFICATIONS


Model	BBWA-3201
Electrical Requirement	220V 50Hz
Power Consumption	500W
Temperature Range	RT+5~100°C
Temperature Stability	±0.5°C
Display Resolution	0.1°C
Chamber Volume	15.8L
Interior Dimension (WxDxH, mm)	330x300x160
Exterior Dimension (WxDxH, mm)	570x220x260
Timing Range	1~9999min
Packing Dimension	65*31*35
CBM	0.07
Net Weight	10.00
Gross weight	12
Alt Name	Water Bath

## ACCESSORIES FOR PURCHASE

No	Name
1	Intelligent programmable temperature controller
2	RS485 connector
3	Printer



1



1



2



2

## APPLICATIONS

Provided for precise and constant temperature and auxiliary heating in colleges industrial and mining enterprises and scientific research departments.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)