



WATER BATH BBWA-3204

WATER BATH BBWA-3204

CIRCULATING BATH, WATER BATH



Microprocessor controller (with timing function)

- Microprocessor control with timing function.
- Digital display.
- Audible and visible alarm for over temperature.

SPECIFICATIONS


Model	BBWA-3204
Electrical Requirement	220V 50Hz
Power Consumption	800W
Temperature Range	RT+5~100°C
Temperature Stability	±0.5°C
Display Resolution	0.1°C
Chamber Volume	21.6L
Interior Dimension (WxDxH, mm)	450x300x160
Exterior Dimension (WxDxH, mm)	460x280x190 / 610x340x260 / 760x340x280
Timing Range	1~9999min
Packing Dimension	70*42*38
CBM	0.11
Net Weight	10.00
Gross weight	16
Alt Name	Water Bath (all stainless steel)

ACCESSORIES FOR PURCHASE

No	Name
1	Intelligent programmable temperature controller
2	RS485 connector
3	Printer



1



1



2



2

APPLICATIONS

Provided for precise and constant temperature and auxiliary heating in colleges industrial and mining enterprises and scientific research departments.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com