



## WATER BATH BBWA-3206

# WATER BATH BBWA-3206

CIRCULATING BATH, WATER BATH



Microprocessor controller (with timing function)

- Microprocessor control with timing function.
- Digital display.
- Audible and visible alarm for over temperature.

## SPECIFICATIONS

Model	BBWA-3206
Electrical Requirement	220V 50Hz
Power Consumption	750W
Temperature Range	RT+5~99°C
Temperature Stability	±0.3°C
Display Resolution	0.1°C
Chamber Volume	2.1Lx3
Interior Dimension (WxDxH, mm)	162x149x151
Exterior Dimension (WxDxH, mm)	505x245x350
Timing Range	1~9999min
Packing Dimension	55*33*40
CBM	0.07
Net Weight	8.00
Gross weight	11
Alt Name	Water Bath (three holes)

## ACCESSORIES FOR PURCHASE

No	Name
1	Intelligent programmable temperature controller
2	RS485 connector
3	Printer



1



1



2



2

## APPLICATIONS

Provided for precise and constant temperature and auxiliary heating in colleges industrial and mining enterprises and scientific research departments.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)