



PRODUCT CATALOG



LABORATORY HIGH SPEED CENTRIFUGE BCBH-307

LABORATORY HIGH SPEED CENTRIFUGE BCBH-307

HEMATOCRIT CENTRIFUGE



1. Brushless DC motor, no carbon dust pollution, free maintenance.
2. Microprocessor control, small LCD display.
3. Automatic rotor identification function.
4. Electronic lock, over speed and imbalance protection.
5. 10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice.
6. It is with CE & ISO9001 & ISO13485 Certificates.

SPECIFICATIONS

Model	BCBH-307
Max. Speed	12000rpm
Max. RCF	15800xg
Max. Volume	24 capillary tubes
Time Range	0~99min59sec
Noise	≤ 55dBA
Speed Accuracy	±20rpm
Size (W x D x Hmm)	335*270*190mm
Net Weight(Kg)	13KG
Power supply	AC 220V 50HZ 2A
Package	Carton box
Alt Name	Hematocrit centrifuge



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com