



PRODUCT CATALOG



LABORATORY HIGH SPEED CENTRIFUGE BCBH-401

LABORATORY HIGH SPEED CENTRIFUGE BCBH-401

LARGE CAPACITY HIGH SPEED CENTRIFUGE



1. Widely used in the field of biochemistry, biological products, pharmaceutical factory and laboratory.
2. Brushless frequency motor for this model which in great torque, free maintenance, no powder pollution, quick in speed up and down.
3. This model with the digital display which indicates the speed,time
4. Micro computer control, there are 10 kinds of program and 10 kinds of acceleration and deceleration for your choice.
5. There are high speed rotors and low speed rotors for your choice.
6. Electric lid lock, compact design, super speed and imbalance protection.
7. The centrifuge body is made of high-quality steel, safe and reliable.

SPECIFICATIONS

| | |
|---------------------|--|
| Model | BCBH-401 |
| Max. Speed | 12000 rpm |
| Max. RCF | 22650 x g |
| Max. Capacity | 4 x 1000 ml |
| Timer | 1~9h59min |
| RPM/RCF Convert | Yes |
| Noise (dB) | ≤ 60 |
| Acc/Dec | 10 Kinds |
| Speed Accuracy | ±20 r/min |
| Voltage(V/Hz) | AC 220V/110V 50HZ/60HZ |
| Size (L x W x H mm) | 780 x 675 x 830 mm |
| Net Weight(Kg) | 200 KG |
| Certificates | CE, ISO & Calibration report are available |
| Alt Name | Large Capacity High Speed Centrifuge |

ACCESSORIES FOR PURCHASE

| No | Name | Max Speed (rpm) | Max Volume (ml) | Max RCF (xg) |
|----|------------------|-----------------|-----------------|--------------|
| 1 | Angle Rotor | 12000 | 6x300ml | 22650 |
| 2 | Angle Rotor | 8000 | 6x500ml | 11740 |
| 3 | Swing Rotor | 4000 | 4x1000ml | 4060 |
| 4 | Microplate Rotor | 4000 | 4x4x96well | 2940 |
| 5 | Angle Rotor | 12000 | 6x250ml | 23670 |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com