



LABORATORY LOW SPEED CENTRIFUGE BCBL-301

LABORATORY LOW SPEED CENTRIFUGE BCBL-301

TABLETOP LOW SPEED CENTRIFUGE(PORTABLE CENTRIFUGE)



Features:

1. Widely used in laboratory, hospital and blood bank.
2. Brushless DC motor for model, free maintenance, no powder pollution, quick in speed up and down.
3. Range of speed from 0 to 4000rpm, smooth in operation, low noise and small vibration.
4. Micro computer control system, digital display in the RCF,time and speed. There are 10 kinds of acceleration and deceleration for your choice.
5. With lid lock, compact design, super speed and imbalance protection.
6. With over speed and imbalance protection, it's safe and reliable

SPECIFICATIONS

Model	BCBL-301
Max. Speed	4000rpm
Max. RCF	1800xg
Max. Capacity	8x15ml
Time Range	0~99min59sec
Noise	≤ 45dBA
Speed Accuracy	±20rpm
Size (W x D x Hmm)	355x320x260mm
Net Weight(Kg)	7KG
Power supply	AC 220V 50HZ 2A
Package	Carton box
Note	Centrifuge+any one rotor
Alt Name	Tabletop low speed centrifuge

ITEMS INCLUDED

No	Rotor no. (Any one rotor)	Name	Max Speed (rpm)	Max Volume (ml)	Max RCF (g)
1	31701	Angle Rotor	4000	8x15ml	1800
2	31702	Angle Rotor	4000	6x10ml/15ml	1800
3	31703	Angle Rotor	4000	12x5ml	1600
4	31704	Angle Rotor	4000	12x7ml	1720
5	32666	Angle Rotor	4000	8x10ml vacuum tube	1800
6	32666	Angle Rotor	4000	8x5ml vacuum tube adapter	1600
7	32666	Angle Rotor	4000	8x7ml vacuum tube adapter	1800



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com