



LABORATORY LOW SPEED CENTRIFUGE BCBL-307

LABORATORY LOW SPEED CENTRIFUGE BCBL-307

CYTOPREP-1 CYTO CENTRIFUGE



1. Brushless DC motor which is in great torque, free maintenance, no powder pollution, quick in speed up and down.
2. Microprocessor-controlled, digital display in speed, time, RCF in operation, speed raising and reducing quickly, operate simply.
3. Automatically electric lid lock, super speed, over temperature protection and imbalance protection.
4. 3 layers protection of steel structure and get the ideal centrifugation result.
5. 10 kinds of accelerating and decelerating speed for your choice.
6. Tabletop centrifuge with metal body, metal lid to ensure the safety of people and instrument.

SPECIFICATIONS

Model	BCBL-307
Max. Speed	4000rpm
Max. RCF	2170xg
Max. Volume	6x9ml
Time Range	0~99min59sec
RPM/RCF Convert	Yes
Noise	≤ 55dBA
Temperature	Normal
Acc/Dec	10 Kinds
Speed Accuracy	±20rpm
Size (W x D x Hmm)	528x370x280mm
Net Weight(Kg)	26KG
Power supply	AC 220V/110V 50HZ/60HZ
Certificates	CE, ISO & Calibration report are available
Alt Name	Cytoprep-1 Cyto Centrifuge



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com