



LABORATORY LOW SPEED REFRIGERATED CENTRIFUGE BCBLR-102

LABORATORY LOW SPEED REFRIGERATED CENTRIFUGE

BCBLR-102

TABLETOP LOW-SPEED REFRIGERATED CENTRIFUGE



It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.

Digital displays make visualized operations convenient.

Imported compressor unit with Non-CFC refrigerant R134a.

It has self-locking, over-speed detection, over-temperature detection, and imbalance detection that is an automatic alarm device.

It has 9 levels acceleration and deceleration to be setting according to the different requirements.

Silastic sealed ring comfort to GMP, US FDA certificated.

SPECIFICATIONS

Model	BCBLR-102
Max speed	8000r/min
Max RCF	9100xg
Max capacity	4x250ml
Speed accuracy	±30r/min
Time setting	1min~99min
Temperature setting	-5°C ~ +35°C
Temperature accuracy	±1°C
Compressor unit	French Tecumseh Compressor Unit R134a
Noise	<65dB(A)
Power supply	AC110V/AC220V, 50Hz~60Hz
Consumption	1.2kW
Dimension (LxWxH)	460mmx710mmx420mm
Packing size (LxWxH)	580mmx730mmx520mm
Net weight (No rotor)	80kg
Alt Name	Tabletop Low-speed Refrigerated Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Description
1	Swing rotor	Capacity: 4x250 ml Max Speed: 4300 r/min Max RCF: 3582 xg
2	Adapter	250 ml Bucket
3	Adapter	250 ml PE Bottle
4	Adapter	4x9x1.5 ml
5	Adapter	4x7x5 ml
6	Adapter	4x9x5/7 ml (Blood Vacuum Tube)
7	Adapter	4x7x10 ml (Internal Thread)
8	Adapter	4x5x15 ml
9	Adapter	4x4x20 ml
10	Adapter	4x50 ml
11	Adapter	4x100 ml

12	Angle Rotor	Capacity: 8x100 ml Max Speed: 8000 r/min Max RCF: 9100 xg Adaptor: 10 ml, 20 ml
13	Swing rotor	Capacity: 4x100 ml Max Speed: 5000 r/min Max RCF: 4390 xg
14	Adapter	4x4x1.5 ml 4x25 ml 4x15 ml 4x50 ml 4x4x5/7 ml 4x3x10 ml 4x20 ml 4x100ml
15	Swing rotor	Capacity: 4x18x7 ml Max Speed: 4000 r/min Max RCF: 2810 xg 75 mm length tube only
16	Swing rotor	Capacity: 4x12x7 ml Max Speed: 4000 r/min Max RCF: 2810 xg 75 mm & 100 mm length tube
17	Swing rotor	Max Speed: 4000r/min Max RCF: 2810xg
18	Adapter	4x4x15 ml 1.5 ml, 5 ml, 7 ml, 10 ml
19	Adapter	4x8x15 ml 1.5 ml, 5 ml, 7 ml, 10 ml
20	Adapter	4x6x15 ml 1.5 ml, 5 ml, 7 ml, 10 ml
21	Adapter	4x2x50 ml 2x1.5 ml 5 ml, 7 ml 10 ml, 15 ml, 20 ml
22	Microplate Bucket	Capacity: 2x2x96 Max Speed: 4000 r/min Max RCF: 2200 xg



1



2



3



4



5



6



7



8



9



10



11



12



13



14



15



16



17



18



19



20



21



22



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com