



PRODUCT CATALOG



LABORATORY LOW SPEED REFRIGERATED CENTRIFUGE BCBLR-303

LABORATORY LOW SPEED REFRIGERATED CENTRIFUGE

BCBLR-303

TABLETOP LOW SPEED REFRIGERATED CENTRIFUGE



1. Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running.
2. Imported compressors, fluorine-free, double cycle cooling, cold and hot alternating easily, environment pollution-free and precise in temperature control.
3. Microprocessor control, LCD color display with touch key screen which indicates the speed, time, RCF in operation, speed raising and reducing quickly, operate simply. Running parameters can be edited.
4. 10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice.
5. Automatically electric lid lock, super speed, over temperature protection and imbalance protection. The centrifuge body is made of high quality steel, safe and reliable.

SPECIFICATIONS

Model	BCBLR-303
Max. Speed	5000rpm
Max. RCF	4730xg
Max. Capacity	4x250ml
Time Range	0~99min59sec
Noise	≤ 60dBA
Speed Accuracy	± 20r/min
Size (W x D x Hmm)	650x665x400mm
Net Weight(Kg)	110KG
Power supply	AC 220V 50HZ 15A
Temperature Range	-20~40°C
Package	Wooden box
Alt Name	Tabletop Low Speed Refrigerated Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max speed (rpm)	Max Volume(ml)	Max RCF (g)
1	Swing out rotor	4000	4x250ml (round cup)	2880
2	Swing out rotor	4000	4x250ml (square cup)	2830
3	Swing out rotor	5000	4x50ml	4730
4	Swing out rotor	5000	4x100ml	205
5	Swing out rotor	4000	8x50ml	189
6	Swing out rotor	4000	16x15ml	3200
7	Swing out rotor	4000	24x15ml	189
8	Swing out rotor	4000	48x5ml vacuum tube	2480
9	Swing out rotor	4000	48x7ml vacuum tube	2760
10	Swing out rotor	4000	40x10ml vacuum tube	2880
11	Swing out rotor	4000	32x10ml vacuum tube	2580
12	Swing out rotor	4000	32x7ml vacuum tube adapter	2580
13	Swing out rotor	4000	32x5ml vacuum tube adapter	2580
14	Swing out rotor	4000	24x10ml vacuum tube	2580
15	Swing out rotor	4000	24x7ml vacuum tube adapter	2580

16	Swing out rotor	4000	24x5ml vacuum tube adapter	2580
17	Swing out rotor	4000	16x10ml vacuum tube	2580
18	Swing out rotor	4000	16x7ml vacuum tube adapter	2580
19	Swing out rotor	4000	16x5ml vacuum tube adapter	2580
20	Micro plate rotor	4000	2x3x48 well	2300
21	Micro plate rotor	4000	2x3x96 well	2490

FEATURES

1. Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running.
2. Imported compressors, fluorine-free, double cycle cooling, cold and hot alternating easily, environment pollution-free and precise in temperature control.
3. Microprocessor control, LCD color display with touch key screen which indicates the speed, time, RCF in operation, speed raising and reducing quickly, operate simply. Running parameters can be edited.
4. 10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice.
5. Automatically electric lid lock, super speed, over temperature protection and imbalance protection. The centrifuge body is made of high quality steel, safe and reliable.
6. Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction.
7. There are several kinds of rotors for your choice, adapters can be designed by experimental requirements.
8. 3 layers protection of steel structure and get the ideal centrifugation result.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
 Email: info@biolabscientific.com | Website: www.biolabscientific.com