



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-103

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

BCBT-103

WEAK ACID AND ALKALI CHEMICALS STORAGE CABINET

The Weak Acid and Alkali Chemicals Storage Cabinet is a specialized storage solution designed to safely house chemicals with low corrosivity, such as weak acids and weak alkalis.



SPECIFICATIONS

Model	BCBT-103
Capacity	22 gallon(83L)
Material	Double-layer anti-fire cold-rolled steel,PP tray
Adjustable Shelf	3pcs
Door	Single/Manual
External Size(W*D*H)mm	600*460*1650
Package Size(W*D*H)mm	680*540*1790
Gross Weight	93kg
Alt Name	Weak Acid and Alkali Chemicals Storage Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com