



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-106

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

BCBT-106

WEAK ACID AND ALKALI CHEMICALS STORAGE CABINET

The Weak Acid and Alkali Chemicals Storage Cabinet is a specialized storage solution designed to safely house chemicals with low corrosivity, such as weak acids and weak alkalis.



SPECIFICATIONS

Model	BCBT-106
Capacity	60 gallon(227L)
Material	Double-layer anti-fire cold-rolled steel,PP tray
Adjustable Shelf	2pcs
Door	Double/Manual
External Size(W*D*H)mm	860*860*1650
Package Size(W*D*H)mm	940*940*1790
Gross Weight	170kg
Alt Name	Weak Acid and Alkali Chemicals Storage Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com