



STRONG ACID AND ALKALI STORAGE CABINET BCBT-202

STRONG ACID AND ALKALI STORAGE CABINET BCBT-202

STRONG ACID & ALKALI STORAGE CABINET

The Strong Acid & Alkali Storage Cabinet is a specialized safety cabinet designed for the safe storage of corrosive chemicals, including strong acids and alkalis, in laboratory and industrial environments.



1. Padlock design.
2. Adjustable shelves.
3. Anti-corrosive materials for strong acid and alkali.
4. Customizable.

SPECIFICATIONS

Model	BCBT-202
Capacity	12 gallon(45L)
Material	PP
Adjustable Shelf	1pc
Door	Single/Manual
External Size(W*D*H)mm	590*460*890
Package Size(W*D*H)mm	675*545*1025
Gross Weight	30kg
Alt Name	Strong Acid & Alkali Storage Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com