



## STRONG ACID AND ALKALI STORAGE CABINET BCBT-206

# STRONG ACID AND ALKALI STORAGE CABINET BCBT-206

## STRONG ACID & ALKALI STORAGE CABINET

The Strong Acid & Alkali Storage Cabinet is a specialized safety cabinet designed for the safe storage of corrosive chemicals, including strong acids and alkalis, in laboratory and industrial environments.



1. Padlock design.
2. Adjustable shelves.
3. Anti-corrosive materials for strong acid and alkali.
4. Customizable.

## SPECIFICATIONS

|                        |                                      |
|------------------------|--------------------------------------|
| Model                  | BCBT-206                             |
| Capacity               | 45 gallon(170L)                      |
| Material               | PP                                   |
| Adjustable Shelf       | 2pcs                                 |
| Door                   | Double/Manual                        |
| External Size(W*D*H)mm | 1090*460*1650                        |
| Package Size(W*D*H)mm  | 1175*545*1785                        |
| Gross Weight           | 72kg                                 |
| Alt Name               | Strong Acid & Alkali Storage Cabinet |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)