



# STRONG ACID AND ALKALI STORAGE CABINET BCBT-208

# STRONG ACID AND ALKALI STORAGE CABINET BCBT-208

## STRONG ACID & ALKALI STORAGE CABINET

The Strong Acid & Alkali Storage Cabinet is a specialized safety cabinet designed for the safe storage of corrosive chemicals, including strong acids and alkalis, in laboratory and industrial environments.



1. Padlock design.
2. Adjustable shelves.
3. Anti-corrosive materials for strong acid and alkali.
4. Customizable.

## SPECIFICATIONS

Model	BCBT-208
Capacity	90 gallon(340L)
Material	PP
Adjustable Shelf	2pcs
Door	Double/Manual
External Size(W*D*H)mm	1090*860*1650
Package Size(W*D*H)mm	1173*943*1785
Gross Weight	107kg
Alt Name	Strong Acid & Alkali Storage Cabinet



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)