



## SLIDE CABINET BCBT-301-A

# SLIDE CABINET BCBT-301-A

SLIDE CABINET



## SPECIFICATIONS

Model	BCBT-301-A
Base Material	SPCC 1.2 mm cold-roll steel sheets (BaoSteel)
Main Body and Drawers	SPCC 0.8 mm cold-roll steel sheets
Configuration	Base + 12 layers of 72 drawers (6 drawers per layer, total 72 drawers)
Capacity	Accommodates up to 64,800 microscope slides (900 per drawer)
Drawer Locking	Each layer can be individually locked
Sliding Rails	ABS sliding rails ensure smooth sliding and hidden locks prevent unexpected opening
Handles	Chrome-coated ABS
Drawer Interior	Each drawer divided into six sections for convenient filing
Surface Treatment	Oil removal, conditioning, zinc phosphate coating, sandblasting, and powder-spray coating
Slot Manufacturing	Manufactured using metal mold
Labeling Slots	Integrated stamping and forming operation
Locking Mechanism	Each drawer can be individually locked; technology under patent application
Dimensions (WxDxH)	450x478x1620 mm
Alt Name	Slide Cabinet



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)