



SLIDE CABINET BCBT-301

SLIDE CABINET BCBT-301

SLIDE CABINET



SPECIFICATIONS

Model	BCBT-301
Base Material	SPCC 1.2 mm cold-rolled steel
Main Body and Drawer Material	SPCC 0.8 mm steel
Drawers	84 drawers (12 layers, 7 per layer)
Capacity	Holds up to 75,600 microscope slides (900 per drawer)
Sliding Rails	ABS sliding rails ensure smooth movement
Locks	Hidden locks prevent accidental opening
Handles	Chrome-coated ABS
Surface Treatment	Oil removal, conditioning, zinc phosphating, sandblasting, and powder coating
Slot Manufacturing	Specifically manufactured using a metal mold
Labeling Slots	Manufactured through integrated stamping and forming operation
Drawer Locking	Each drawer can be individually locked
Dimensions (WxDxH)	450x478x1620 mm
Alt Name	Slide Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com