



# PRODUCT CATALOG



## LABORATORY LOW SPEED REFRIGERATED CENTRIFUGE BCFHR-305

# LABORATORY LOW SPEED REFRIGERATED CENTRIFUGE

## BCFHR-305

HIGH SPEED REFRIGERATED CENTRIFUGE



1. Widely used in the field of pharmaceutical factory and laboratory.
2. Brushless frequency motor, in great torque, free maintenance, no powder pollution, quick in speed up and down.
3. Digital display which indicates the program, speed, time, RCF, temperature.
4. Micro computer control, there are 10 kinds of program and 10 kinds of acceleration and deceleration for your choice.
5. Several rotors for your choice.
6. Electric lid lock, compact design, super speed and imbalance protection.
7. The centrifuge body is made of high-quality steel, safe and reliable.

## SPECIFICATIONS

Model	BCFHR-305
Max speed	1000 rpm
Max volume	6x1000ml
Timer	1~9h59min
Dimension(HxDxW)	800x740x930mm
Speed accuracy	±20rpm
Temperature range	-20~40°C
Max RCF	15730xg
Noise	≤65dB
Net weight	360KG
Power supply	AC 220V 50HZ/60HZ
Temperature accuracy	±1°C
RPM/RCF Convert	Yes
Acc/Dec	10 Kinds
Certificates	CE, ISO & Calibration report are available
Alt Name	Floor Type Low Speed Refrigerated Centrifuge

## ACCESSORIES FOR PURCHASE

No	Name	Max speed(rpm)	Max Volume(ml)	Max RCF(xg)
1	Angle rotor	8000	6x1000ml	14550
2	Angle rotor	9000	6x500ml	14860
3	Angle rotor	10000	6x300ml	15730



## **Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)