



LABORATORY LOW SPEED CENTRIFUGE BCFL-402

LABORATORY LOW SPEED CENTRIFUGE BCFL-402

FLOOR STANDING LARGE CAPACITY CENTRIFUGE



1. Widely used in the field of blood bank, pharmaceutical factory and laboratory.
2. Brushless frequency motor, in great torque, free maintenance, no powder pollution, quick in speed up and down.
3. Digital display which indicates the program, speed, time, RCF, temperature.
4. Micro computer control, there are 10 kinds of program and 10 kinds of acceleration and deceleration for your choice.
5. Electric lid lock, compact design, super speed and imbalance protection.
6. The centrifuge body is made of high-quality steel, safe and reliable.

SPECIFICATIONS

| | |
|----------------------|--|
| Model | BCFL-402 |
| Max. Speed | 6000 rpm |
| Max. RCF | 6600 x g |
| Max. Capacity | 4 x 1000 ml |
| Timer | 1~9h59min |
| RPM/RCF Convert | Yes |
| Noise (dB) | ≤ 65 |
| Temperature Range | / |
| Acc/Dec | 10 Kinds |
| Speed Accuracy | ±20 r/min |
| Temperature Accuracy | / |
| Voltage(V/Hz) | AC 220V/110V 50HZ/60HZ |
| Size (L x W x H mm) | 800 x 645 x 830 mm |
| Net Weight(Kg) | 170 KG |
| Certificates | CE, ISO & Calibration report are available |
| Alt Name | Floor Standing Large Capacity Centrifuge |

ACCESSORIES FOR PURCHASE

| No | Name | Max Speed (rpm) | Max Volume (ml) | Max RCF (xg) |
|----|-------------|-----------------|-------------------------|--------------|
| 1 | Swing Rotor | 4000 | 4x1000ml (round bottle) | 4060 |
| 2 | Angle Rotor | 6000 | 4x300ml | 5390 |
| 3 | Angle Rotor | 6000 | 6x300ml | 5660 |
| 4 | Angle Rotor | 6000 | 6x500ml | 6600 |
| 5 | Microplate | 4000 | 4x4x96well | 2840 |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com