



COMPOUNDING HOOD BCHL-403

COMPOUNDING HOOD BCHL-403

COMPOUNDING HOOD




- There are two types:
- Positive pressure in work area only protects sample
 - Negative pressure in work area protects operator and environment.

SPECIFICATIONS


Model	BCHL-403
External Size (W*D*H)mm	700*620*1150mm
Internal Size (W*D*H)	640*550*700mm
Display	LED Display
Airflow Velocity	0.3~0.5m/s
Material Main body	Cold-steel with anti-bacterial powder coating
Material Work table	304 stainless steel
Pre-Filter	Polyester fiber, washable
Filter	HEPA Filter
Noise	≤60dB
Front Window	Manual, acrylic material
Max Opening	660mm
LED Lamp	4W*2
UV Lamp	15W*1 Emission of 253.7 nanometers
Consumption	300W
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz;
Standard Accessory	HEPA filter, LED lamp*2, UV lamp*2.
Optional Accessory	Base stand
Gross Weight	100kg
Package Size (W*D*H)mm	840*760*1400
Alt Name	Compounding Hood

FEATURES

There are two types:



- Positive pressure in work area only protects sample
- Negative pressure in work area protects operator and environment.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com