



CHROMATOGRAPHY BCHR-104

CHROMATOGRAPHY BCHR-104

ION CHROMATOGRAPHY

Chromatography is a technique that enables the separation, identification, and purification of the components of a mixture for qualitative and quantitative analysis. Our extensive range offers a variety of products like Gas, Ion and Portable Ion chromatography products to meet all separation needs, including improved resolution, enhanced sensitivity, faster analysis and consistent performance.



Built-in eluent generator, free from configuring eluent, with gradient elution available.
Modular manufacturing process to maintain excellent systemic stability.
Built-in low-pressure degassing technology to eliminate bubble interference for more stability.
Optional intelligent automatic injection system for large sample volumes, which features automatic dilution to save labor and time.
Work across a variety of detectors, to expand the scope of applications of ion chromatography.

SPECIFICATIONS

Model	BCHR-104
Ion Chromatographic Pump	
- Maximum Pressure	35 MPa (PEEK)
- Type	High-pressure and low-pulse two-piston tandem advection pump
- Pressure Display Accuracy	≤ 0.1 MPa
- Flow Range	0.001 ~ 9.999 mL/min
- Pressure Pulse	$\leq 0.5\%$
- Flow Stability	(0.2-0.5) mL/min $\leq 3\%$; (0.5-1.0) mL/min $\leq 2\%$; >1.0 mL/min $\leq 2\%$
Conduction Detection System	
- Cell Volume	$\leq 0.8\mu\text{L}$
- Detection Mode	Bipolar conductivity detection
- Detection Range	0~50000 $\mu\text{S/cm}$
- Detection Resolution	$\leq 0.0020\text{nS/cm}$
- Output Voltage	-6000~+6000 mV (adjustable)
- Electronic Noise	0.02 nS
- Baseline Noise	$\leq 0.05\%$ FS
- Baseline Drift	$\leq 3\%$ FS
- Operating Temperature Range	Room temperature +5°C~60°C (41~140°F)
- Controlling Temperature Accuracy	$\pm 0.01^\circ\text{C}$
- Linear Range	$\geq 10^3$
- Instrument Linearity	≥ 0.999
- Quantitative Repeatability	$\leq 1.0\%$
- Qualitative Repeatability	$\leq 1.0\%$
- Minimum Detectable Concentration	$\text{Cl}^- \leq 0.0005 \mu\text{g/mL}$; $\text{Li}^+ \leq 0.001 \mu\text{g/mL}$; $\text{BrO}_3^- \leq 0.001 \mu\text{g/mL}$
Flow System	
- Six-way Valve	PEEK material, pressure 5000 psi; Independent automatic collecting and flow function.
Built-in Eluent Generator	

- Eluent Types	KOH/MSA
- Eluent Concentration Range	0.1-120 mM
- Concentration Increment	0.1 mM
- Flow Rate Range	0.1-5.0 mL/min
- Minimum Pressure	5 MPa
Suppressor	
- Dead Volume	<50 µL
Other Specifications	
- Dimension (LxWxH)	350x470x650 mm
- Net Weight	34 kg
- Gross Weight	40 kg
- Power	150 W
Alt Name	Ion Chromatography

APPLICATIONS

Food Testing, Chemical Industry, Beverage Testing, Drug testing, Forensic Science, Pharmaceutical, Molecular Biology, Medical, Research, Laboratory.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com