



## HIGH TEMPERATURE CIRCULATOR BCHT-104

# HIGH TEMPERATURE CIRCULATOR BCHT-104

## HIGH TEMPERATURE CIRCULATOR

We offer you robust circulators for high temperature applications equipped with a circulating pump to provide excellent temperature uniformity and safe operation. Easy to maintain with LCD display, these circulators are user friendly and give reliable results.



Microprocessor temperature controller  
Easy maintenance  
LCD display with touch key panel for temperature setting  
Over temperature alarm system

## SPECIFICATIONS

Model	BCHT-104
Capacity	30 L
Temperature Range	RT+5°C to 100°C
Temperature Accuracy	±0.05°C
Pump Flow Rate	6 L/min
Top Opening	310x280 mm
Inner Dimension (WxDxH)	400x330x230 mm
Outer Dimension (WxDxH)	550x490x600 mm
Weight	21.62 kg
Alt Name	High Temperature Circulator

## APPLICATIONS

Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)