



MINI CENTRIFUGE BCMI-101

MINI CENTRIFUGE BCMI-101

MINI CENTRIFUGE



SPECIFICATIONS

Model	BCMI-101
Max speed	4000r/min
Max RCF	900xg
Max capacity	6x1.5ml/0.5ml/0.2ml
Noise	<48dB(A)
Power supply	AC110V~AC220V
Consumption	30W
Dimension (LxWxH)	189mmx161mmx122mm
Packing size (LxWxH)	330mmx250mmx210mm
Net weight (No rotor)	1kg
Alt Name	Mini Centrifuge



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com