



MINI CENTRIFUGE BCMI-102

MINI CENTRIFUGE BCMI-102

MINI CENTRIFUGE



SPECIFICATIONS

| | |
|-----------------------|-------------------|
| Model | BCMI-102 |
| Max speed | 6000r/min |
| Max RCF | 2000xg |
| Max capacity | 8x0.2ml |
| Noise | <48dB(A) |
| Power supply | AC110V~AC220V |
| Consumption | 30W |
| Dimension (LxWxH) | 189mmx161mmx122mm |
| Packing size (LxWxH) | 330mmx250mmx210mm |
| Net weight (No rotor) | 1kg |
| Alt Name | Mini Centrifuge |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com